

## A Novel RT-LAMP Assay for Rapid and Simple Detection of Classical Swine Fever Virus\*

Lei CHEN, Xue-zheng FAN, Qin WANG\*\*, Lu XU, Qi-zu ZHAO,  
Yuan-chen ZHOU, Jun LIU, Bo TANG and Xing-qi ZOU  
(China Institute of Veterinary Drug Control, Beijing, 100081, China)

**Abstract:** A simple and rapid assay for the detection of Classical swine fever virus (CSFV) was established using reverse transcription loop-mediated isothermal amplification (RT-LAMP). This study describes the amplification of the genomic RNA of CSFV under isothermal conditions (63°C) within one hour, using a set of six primers (two outer primers, two inner primers and two loop primers). This RT-LAMP assay showed 100-fold higher sensitivity than the standard RT-PCR method and identified eighteen additional positive cases that were negative when tested by RT-PCR. This RT-LAMP was able to detect all the 13 strains of CSFV but not the BVDV, PRRSV, SIV, PRV-PCV, thus showed a good specificity. Products amplified by RT-LAMP can be visualized by agarose gel electrophoresis and in addition, either as a white precipitate at the bottom of the tube after a pulse spin or as a color change when dyed with SYBR Green I which are visible to the naked eye. Because RT-LAMP is low-cost and produces rapid results, it has the potential to be an excellent tool for CSFV surveillance in the field, especially in developing countries.

**Key words:** Classical swine fever virus (CSFV); Reverse transcription loop-mediated isothermal amplification (RT-LAMP); Rapid detection

Classical swine fever virus (CSFV), in the *Pestivirus* genus and the *Flaviviridae* family, is a small, enveloped, positive and single-stranded RNA virus. CSFV is the causative agent of Classical swine fever (CSF). CSF, also known as ‘hog cholera’, is a highly

contagious febrile disease worldwide. Outbreaks of CSF cause heavy losses in pig production and severely hamper the international trade in livestock. For these reasons, CSF is listed as a notifiable disease by the Office International des Epizooties and the European Union<sup>[1, 9]</sup>. Thus early diagnosis of the disease is important to effectively control the spread of the CSFV.

Virus isolation is the ‘gold standard’ for CSFV diagnosis. However, it cannot meet the rapid detection requirements because it is labor intensive and time

Received: 2009-03-30, Accepted: 2009-10-27

\* Foundation items: The National Science and Technology supporting plan of the Eleventh Five-year (2006BAD06A18 and 2006BAD06A03); Beijing Natural Science Foundation (5072041).

\*\* Corresponding author.

Phone: +86-10-62103680, Fax: +86-10-62103674,

E-mail: wangqin@ivdc.gov.cn

consuming. Direct fluorescent antibody test of frozen sections is rapid and reliable, but needs to be performed by well-trained technicians. Antigen-capture enzyme-linked immunosorbent assays are suitable for early diagnosis but are less sensitive. Recent molecular techniques based on genomic sequence detection such as reverse transcription PCR (RT-PCR), reverse transcription nested PCR (RT-nPCR) and real-time PCR are sensitive and rapid, but these methods use thermocycling equipment and so are impractical for field applications and for use in less developed laboratories [5].

Recently, a novel, rapid and sensitive technique named loop-mediated isothermal amplification (LAMP) was developed by Notomi *et al* [8]. It can amplify a target gene (with a DNA or an RNA template) in a very short time with high sensitivity, specificity and efficiency under isothermal conditions [6, 7, 10, 13]. For these reasons, LAMP has been widely used for virus diagnosis [2, 3, 11, 14]. In this study, a one-step, single-tube RT-LAMP assay was developed for the rapid detection of CSFV.

## MATERIALS AND METHODS

### Viruses isolates

The reference strains CSFV SM, HCLV and Thiverval strains were provided by the China Veterinary Culture Collection Center. The other CSFV strains BJ1/08, BJ2/08, JL1/94(blood), GXBH1/98, HeBZJK1/99, BJCY1/96, USA-331, JL1/94(cell culture), GDGZ1/95 and FJFQ1/98 were isolated and maintained by the National Classical Swine Fever Reference Laboratory, China. Bovine viral diarrhoea virus (BVDV), Porcine circovirus (PCV), Porcine reproductive and Respiratory syndrome virus (PRRSV), Swine influenza virus (SIV), Porcine parvovirus (PPV) and Pseudorabies

virus (PRV) were also obtained from the China Veterinary Culture Collection Center.

### Clinical samples

A total of 157 acute-phase samples, including blood and tonsil, collected from dead pigs with a clinical diagnosis of Classical swine fever and admitted to the National Classical Swine Fever Reference Laboratory, China, were used to evaluate the RT-LAMP method by comparison with RT-PCR. All the samples were transported to the laboratory under strict cold chain conditions and stored at -80°C until further investigation.

### RNA and DNA extractions

Total RNA was extracted with Trizol reagent (Invitrogen, USA) according to the manufacturer's instructions. After elution with 20 µL nuclease-free double-distilled H<sub>2</sub>O, RNA samples were stored at -70°C until required. Total DNA was extracted using the TIANamp Genomic DNA Kit (Tiangen Biotech, Beijing, China) according to the manufacturer's instructions. After extraction, DNA was eluted with 20 µL elution buffer and stored at -20°C for analysis.

### Design of CSFV-specific RT-LAMP primers

To find conserved regions of the CSFV genome, 28 genome sequences of CSFV (GenBank accession numbers: A16790, AF091661, AF092448, AF326963, AF333000, AF352565, AF407339, AY072924, AY367767, AY554397, AY568569, AY578687, AY578688, AY646427, AY663656, AY775178, AY805221, D49532, D49533, DQ127910, J04358, L49347, M31768, NC\_002657, U45478, U90951, X87939, X96550) were downloaded from GenBank and aligned using the Vector NTI 10 program. Several sets of six primers (two outer, two inner and two loop primers) were designed with the Primer Explorer V4 software (<https://primerexplorer.jp>) according to these potential target regions and one




Fig.1. Design of the RT-LAMP for CSFV detection. A: Positions of the oligonucleotide primers and restriction enzyme. The primer sites are indicated by the arrows in the sequence and the restriction enzyme site of *Hinf* I is highlighted in dark grey. B: Sequences of selected primers.

set of primers was selected based on the criteria described by Notomi *et al* (Fig.1)<sup>[8]</sup>.

**RT-PCR**

In order to test the sensitivity of the RT-LAMP assay, RT-PCR was performed using a CSFV RT-PCR kit (Anheal Laboratories, Beijing, China). The same RNA templates were used for both RT-LAMP and RT-PCR. Positive and negative controls were used for each reaction and all reactions were carefully monitored for cross-contamination.

**RT-LAMP assay**

The RT-LAMP reaction was carried out in a total volume of 25  $\mu$ L, with mixture of 1 $\times$ Thermo buffer (New England Biolabs, USA), 0.8 mol/L Betaine (Promega, USA), 1.4 mmol/L dNTP mix (Dingguo, China), 4mmol/L MgSO<sub>4</sub> (Promega, USA), 8 U *Bst* DNA polymerase (New England Biolabs, USA), 2 U

AMV reverse transcriptase (Promega, USA), 5 pmol/L F3 and 5 pmol/L B3, 50 pmol/L FIP and 50 pmol/L BIP, 25 pmol/L LF and 25 pmol/L LB and 2  $\mu$ L of template RNA. The mixture was incubated at 63  $^{\circ}$ C for 60 min, and then heated at 80  $^{\circ}$ C for 10 min to stop the reaction. RT-LAMP assay products (10  $\mu$ L) were electrophoresed on a 2% agarose gel in Tris-borate buffer. The result was visualized with a UV transilluminator at 302 nm. The amplification should be seen as a ladder-like pattern on the gel due to the formation of a mixture of stem-loop DNAs with various stem lengths and cauliflower-like structures, with multiple loops induced by annealing between alternately inverted repeats of the target sequence in the same strand.

To facilitate field application of the RT-LAMP assay, we also monitored the reaction visually. A positive

amplification was seen as either the precipitation of the white by-product magnesium pyrophosphate in the tube bottom after a pulse spin [13] or by a color change following the addition of 1  $\mu$ L of SYBR Green I dye (Invitrogen, USA) to the tube.

### Specificity tests

To confirm the specificity of the RT-LAMP assay, cross-reactivity tests were conducted using different sources of RNA or DNA as templates. Restriction endonuclease digestion was used simultaneously to further confirm the specificity.

## RESULTS

### Comparison of the sensitivity of RT-LAMP with that of RT-PCR

The sensitivity of the CSFV RT-LAMP assay was determined using 10-fold serial dilutions of RNA extracted from 0.2 mL of  $10^{3.84}$ TCID<sub>50</sub> blood samples




Fig. 2. Comparative sensitivity of RT-LAMP and RT-PCR. RT-LAMP products(A) and RT-PCR products(B). M, 2000-bp DNA marker; 1-8, RT-LAMP products or RT-PCR products of dilution RNAs extracted from the 200  $\mu$ L blood stock of CSFV-Shimen ( $10^{3.84}$  TCID<sub>50</sub>/mL). 1 lane, Template concentration  $10^{-1}$ ; 2 lane,  $10^{-2}$ ; 3 lane,  $10^{-3}$ ; 4 lane,  $10^{-4}$ ; 5 lane,  $10^{-5}$ ; 6 lane,  $10^{-6}$ ; 7 lane,  $10^{-7}$ ; 8 lane,  $10^{-8}$ .

containing the SM virus. The results showed that the accelerated RT-LAMP assay was 100-fold more sensitive compared with conventional RT-PCR, the detection limit of RT-LAMP being a  $10^{-6}$  dilution of template RNA (Fig.2).

### Specificity of CSFV RT-LAMP assay

The RT-LAMP assay detected 13 different isolates of CSFV, each showing the characteristic ladder-like pattern in the gel, the formation of a white precipitate in the reaction tube and a color change after the addition of SYBR Green I dye. Also, the cross-reactivity tests for all non-Classical swine fever viruses were negative (Fig.3). Finally, the amplified products were confirmed by digestion with restriction enzyme *Hinf* I (Fig.4).

### Evaluation of CSFV RT-LAMP assay

A total of 157 samples were used in this study for comparative evaluation between RT-LAMP and RT-PCR. The detection rates of RT-LAMP and RT-PCR were 26.11% and 14.65% respectively. None of the RT-PCR-positive samples were missed by RT-LAMP. Critically, RT-LAMP identified eighteen additional




Fig. 3. Specificity of the CSFV RT-LAMP method. M, 2000-bp DNA marker; 1, RT-LAMP product of CSFV-Shimen; 2-7 lane, RT-LAMP product of Bovine viral diarrhea virus (BVDV), Respiratory syndrome virus (PRRSV), Swine influenza virus (SIV), Porcine parvovirus (PPV), Pseudorabies virus (PRV) and Porcine circovirus (PCV), respectively.




Fig. 4. Restriction analysis of the CSFV RT-LAMP products on 2% agarose gel. M, 2000-bp DNA marker; 1, RT-LAMP product of CSFV-Shimen; 2, RT-LAMP product of CSFV-Shimen digested with *Hinf* I.

Table 1. Comparative evaluation of RT-LAMP assay with RT-PCR for 157 samples with CSF clinical signs.

Specimen type	Tested samples	Positive samples tested by	
		RT-LAMP	RT-PCR <sup>a</sup>
Tonsil	97	36	23
Blood	60	5	0

<sup>a</sup> All positive samples detected by RT-PCR were confirmed by sequencing.

positive cases that were found to be negative by RT-PCR suggesting a higher sensitivity for this assay (Table 1).

## DISCUSSION

Classical swine fever is an important disease in pigs, which causes significant economical losses in the pig industry worldwide and though some countries had supposedly eradicated CSF, it has re-appeared in some areas and led to the slaughter of millions of pigs<sup>[4]</sup>. More importantly, the clinical prevalence of apparent infection with CSFV posed a greater risk to the pig industry. However, current methods for the detection of CSFV such as real-time PCR and RT-nest PCR are either time consuming or require special equipment.

The RT-LAMP assay is a useful diagnostic technique

with many advantages. Firstly, the RT-LAMP reaction takes place under isothermal conditions so no heating equipment is required: a simple water bath or heat pack can accomplish the amplification. Secondly, the assay is simple. The whole reaction is carried out in a single tube by mixing several reagents together at 63°C for 1 h. Thirdly, it is efficient. Large amounts of the products and by-products that RT-LAMP yields can be identified visually. Fourthly, the RT-LAMP is a one-step assay in which the reverse transcription and LAMP reaction are combined by adding AMV reverse transcriptase to the same tube, which further simplifies the assay when compared with the conventional RT-PCR method. Therefore, in the field of disease diagnosis the RT-LAMP assay has the potential to replace current nucleic acid amplification techniques<sup>[2, 12]</sup>.

In this study an RT-LAMP assay was established for CSFV detection. A set of six primers was designed against the conserved coding regions in CSFV and strain SM was used for optimization of the assay. RT-LAMP was able to detect 13 strains of CSFV and showed good specificity. RT-LAMP was able to detect a 10<sup>-6</sup> dilution of RNA extracted from blood with 10<sup>3.84</sup>TCID<sub>50</sub> CSFV, which was 100-fold more sensitive than RT-PCR. For clinical diagnosis of CSFV, the data presented in this study shows that the RT-LAMP assay is more sensitive than RT-PCR as it identified eighteen additional cases that were found to be negative by RT-PCR.

In conclusion, the RT-LAMP assay provides a simple, rapid, sensitive and accurate method for detection of CSFV in clinical samples. The assay does not require any sophisticated equipment and in many situations is a far more suitable technique for the detection of CSFV in the field.

## Reference

1. **Anonymous.** Council Directive 82/894/EEC of 21 December 1982 on the notification of animal diseases within the Community. **Official Journal L**, 378: 58-62.
2. **Dukes J P, King D P, Alexandersen S.** 2006. Novel reverse transcription loop-mediated isothermal amplification for rapid detection of foot-and-mouth disease virus. **Arch. Virol**, 151: 1093-1106.
3. **Fukuda S, Takao S, Kuwayama M, et al.** 2006. Rapid detection of norovirus from fecal specimens by real-time reverse transcription-loop-mediated isothermal amplification assay. **J Clin Microbiol.** 44: 1376-1381.
4. **Greiser-Wilke I, Blome S, Moenning V.** 2007. Diagnostic methods for detection of Classical swine fever virus: *Status quo* and new developments. **Vaccine**, 25: 5524-5530.
5. **Marie F L P, Alain M, Philippe V.** 2006. Classical swine fever and other Pestiviruses. In: **Diseases of Swine, the 9<sup>th</sup> edition** (Straw B E, Zimmerman J, D'Allaire S, et al.). Blackwell Publishing, p316-317.
6. **Mori Y, Nagamine K, Tomita N, et al.** 2001. Detection of loop-mediated isothermal amplification reaction by turbidity derived from magnesium pyrophosphate formation. **Biochem Biophys Res Commun**, 289: 150-154.
7. **Nagamine K, Hase T, Notomi T.** 2002. Accelerated reaction by loop-mediated isothermal amplification using loop primers. **Mol Cell Probes**, 16: 223-229.
8. **Notomi T, Okayama H, Masubuchi H, et al.** 2000. Loop-mediated isothermal amplification of DNA. **Nucleic Acids Res**, 28, e63.
9. **Office International des Epizooties,** 2006. Diseases notifiable to the OIE; 2006, [http://www.oie.int/Eng/maladies/en\\_classification.htm](http://www.oie.int/Eng/maladies/en_classification.htm).
10. **Parida M M, Posadas G, Inoue S, et al.** 2004. Real-time reverse transcription loop-mediated isothermal amplification for rapid detection of West Nile virus. **J Clin Microbiol**, 42: 257-263.
11. **Parida M M, Horioko K, Ishida H, et al.** 2005. Rapid detection and differentiation of dengue virus serotypes by real-time reverse transcription loop-mediated isothermal amplification assay. **J Clin Microbiol**, 43: 2895-2903.
12. **Parida M M, Santhosh S R, Dash P K, et al.** 2006. Development and evaluation of reverse transcription loop-mediated isothermal amplification assay for rapid and real-time detection of Japanese encephalitis virus. **J Clin Microbiol**, 44 (11): 4172-4178.
13. **Savan R, Kono T, Itami T, et al.** 2005. Loop-mediated isothermal amplification: an emerging technology for detection of fish and shellfish pathogens. **J. Fish. Dis.** 28: 573-581.
14. **Toriniwa H, Komiya T.** 2006. Rapid detection and quantification of Japanese encephalitis virus by real-time reverse transcription loop-mediated isothermal amplification. **Microbiol. Immunol.** 50: 379-387.